

Formidabel First

Side protection Unique



Side protection Unique - protection and mobility

As much protection as necessary, as much freedom as possible. Unlike conventional side rails, our Unique side protection offers variable solutions without tools, if required, in just a few steps. Depending on the care situation, the elements can be attached to the bed in different combinations. Side panels when no side guard is needed, side guard elements for residents in need of protection. The height of the element is infinitely adjustable by telescopic extension, so that it protects against falling out when attached to the head or provides support when getting up.



Made in Germany



Formidabel First

Side protection Unique



Formidabel Technical Data

Parameters	Specifications
Overall Width	103cm
Overall Length	210cm
Extended Length	215cm or 225cm
Lying Surface	90 x 200 cm
Siderails	Unique 1/2 + 1/2
Castor Dimension	Ø50mm (Options: Ø75mm, Ø100mm)
Height Adjustment	23cm to 80cm (Caster Ø50mm) 26cm to 83cm (Caster Ø75mm or Ø100mm)
Trendelenburg Adjustment	12°
Backrest Adjustment	70° +/- 5°
Leg Rest Adjustment	9°
Brake System	Castors Ø50mm: 2 rear with single brake + 2 front with single brake (Option for castors Ø75mm or Ø100mm: 4 castors with a central brake)
Safe Working Load	225kg
Mattress Thickness	≤ 20kg
Handset	Scenario
Power Supply	100 to 240 Vac
Backup Battery	3 x 9V Block battery with pole protection